

LEMIS 2000

Net 8 Easy Config to an existing network if not connecting to database

Check another machine on the network

1. *Modify*
2. *(choose from the existing services) check another machine on network – click next*
3. *TCP/IP – click next*
4. *(host name) which is the name of the server or the ip address of the server*
5. *port number is 1521 – click next*
6. *(database sid) – click next*
7. *Test service*
8. *Username: , Password: , for the database - click test*
9. *Successful – click next*
10. *Finish*

Net 8 Easy Config on new workstation

1. *Add new service*
2. *(new service name, usually lemis) check another machine on the network. This can be what you want to name it, this will be what you type in at the 3rd box while logging in to LEMIS – click next*
3. *TCP/IP – click next*
4. *(host name) which is the name of the server or the ip address of the server*
5. *port number: 1521 – click next*
6. *(database sid) – click next*
7. *Test service – click next*
8. *Username: , Password: , for the database – click test*
9. *Successful – click next*
10. *Finish*

****The service name is what you type in the 3rd box (database) while logging into LEMIS 2000. This service name is an alias that points to the SID.*

****The SID can be found on the server listed as OracleServices____ in the services section.*